



Buret, Acrylic, 50 mL

Kat. číslo produktu:	1405941
CZK Cena bez DPH:	Kontaktujte nás
Datum odeslání neuvedeno	

Durable, crystal-clear acrylic body is individually calibrated to meet the accuracy requirements of ASTM E287. Leak-proof stopcock with self-lubricating PTFE TFE plug. Excellent for general lab or field use in schools, industrial, municipal or food applications. Unaffected by dilute mineral acids and bases (except hydrofluoric acid and ammonium hydroxide). Can contain liquids within a pH range of 1 - 14 and up to 5% wt/wt of acids or 30% bases can be used. Not designed for exposure to alcohol or organic solvents.

Technické údaje

Kapacita:	50 mL
Materiál skříně:	Acrylic
Speciální funkce:	Buret
Stupnice:	0.1 mL