



Test kit, carbon dioxide, model CA-DT

Kat. číslo produktu:

2064100

Nebezpečné

CZK Cena bez DPH:

Kontaktujte nás

Datum odeslání neuvedeno



Utilizing the Hach Digital Titrator, the Model CA-DT kit is used to determine carbon dioxide concentrations. Kit includes Digital Titrator, reagents and labware, manual and carrying case. Range: 10.1-1,000 mg/L. Contains approximately 50 tests.

Hach Digital Titrator offers higher levels of precision and accuracy than drop count procedures, with reproducibility comparable to a buret titration

Rugged titration device is ideal for field use, with no glass burets required

Technické údaje

Metoda:	Digital Titrator/ Sodium hydroxide Phenolphthalein
Parametr:	Carbon Dioxide - As CO ₂
Platforma:	Digital Titrator
Počet testů:	50
Rozsah měření:	10 - 1,000 mg/L CO ₂
Smallest increment:	0.1, 2.0
Styl kufříku:	F
Typ testovací soupravy:	CA-DT
Vysvětlivka:	*mg/L unless otherwise noted, ppb = µg/L, ppm = mg/L.; gpg = grains per gallon; 1 gpg = 17.1 mg/L or 17.1 ppm. 1EPA-accepted for wastewater. Check with your regulatory agency.