



Crucible, Gooch, 25 mL Capacity

Kat. číslo produktu:

2098640

CZK Cena bez DPH:

Kontaktujte nás

Datum odeslání neuvedeno

Perforated bottom, glazed inside and out except for outside bottom surface. Filter for Gooch Crucible available and sold separately.

25 mL Gooch Crucible is recommended for use in Standard Methods 2540 C and 2540 D to determine Total Suspended Solids.

Technické údaje

Kapacita: 25 mL

Obsah balení: 1 /pk