



Be Right™



Bottle, Sampling, Polycarbonate, 500 mL, 24/pk

Kat. číslo produktu:	2324579
CZK Cena bez DPH:	Kontaktujte nás
Datum odeslání neuvedeno	

Transparent polycarbonate bottle with cap features wide mouth (3.8 cm/1.48" diameter) opening for easy filling and square design for maximum space efficiency. Autoclavable with molded-in graduations (do not engage threads for autoclave cycle.)

Technické údaje

Body Material:	Polycarbonate
Kapacita:	500 mL
Obsah balení:	24 /pk